



# MEDICAL WIRE BASKET INSTALLATION

ICWUSA.com, Inc.

If you have any questions please call Toll Free 1-800-558-4435

When mounting to studs, ICW recommends using only wood or metal studs. This is the basic method to mount the arm on a stud. Check building code regulations in your area for additional requirements.

## MOUNT TRACK TO WALL


### MD TRACK

1. FOR WOOD STUD
2. **Center mounting option**
3. Locate the wall stud and mark the top center hole of the wall track.  
Drill a 11/64 pilot hole, insert a #14 screw (optional) and tighten slightly. Using a level, mark the lower holes in the center of the track, drill pilot holes and then insert the #14 screws. Use all screws.
4. **Right or left mounting option**
  - Locate the wall stud. Align either right or left side of track with stud and mark one hole in the upper corner for track placement. Drill 11/64" pilot hole. Insert #14 screw and tighten slightly.
  - Using a level, mark the remaining holes on the same side, drill pilot holes into stud. Insert #14 screws and, with the unit level, tighten screws.
  - Insert provided Snap Toggles into outer corner holes on opposite side. See below.
5. FOR METAL STUD
6. #14 **wood** screws are available upon request.  
If installing into **metal studs** you will need to replace these with proper metal stud fasteners or provided snap toggles. All hardware used to secure track mount must have a proper countersunk head.
7. Attach beauty cap.

## IF ATTACHING DIRECTLY TO WALL

### ATTACH TO WALL

1. Recommend using a level.
2. Drill two 2 1/2" screws through bracket holes and into the wall.



NOTICE: ICWUSA.COM, Inc. provides the appropriate hardware fasteners to match the majority if not all standard monitors and devices. It is the responsibility of the installer / integrator to ensure that the screws and fasteners have the adequate thread engagement. The rule of thumb is that the thread depth is greater than the diameter of the hole or 4 full turns of 360 degrees of thread engagement. It is also the responsibility of the installer / integrator to ensure that they do not damage the monitor or device by using a screw or fastener that is too long for hole and therefore damaging the unit by inserting too far.

## IF ATTACHING TO WALL TRACK

### ATTACH TO WALL TRACK

After Wall Track has been installed:

1. Set one (1/4-20 x 7/8) screw (Item 2) through the mounting bracket's top hole in advance.  
Do not tighten yet.
2. Insert track adaptor into track and slide down to the approximate desired position.
3. While holding the screw in the mounting bracket, position wire basket into slots on back side of mounting bracket
4. Align screw with threaded hole in Track Adaptor and tighten several turns by hand. Do not fully tighten yet. Insert bottom screw (Item 2) into Mounting Bracket
5. Align bottom screw with bottom thread hole and tighten by hand. Do not fully tighten yet.
6. Adjust position of basket and slots as necessary to get desired basket placement.
7. Finish tightening with supplied 5/32" hex key (Item 6).




FIGURE 1




FIGURE 2





FIGURE 3



## CLEANING OF THE ICW MOUNTING DEVICE

When cleaning your ICW mounting unit please follow the points below. Proper cleaning of the mounting device will help assure proper function and extend the life of the unit.

1. Any general cleaner is acceptable except degreasers are not recommended near the bearing entrances.
2. Spray cleaner or water on a soft cloth or paper towel and then wipe down the mounts.
3. Do not spray directly at the bearing seams (located at all swivel points) as this could reduce the life of the bearings.