


ICW Custom Medical Device Mount



UL180-BV-W3-A1-DB2NK
The Ultra 180 arm configured with a Drop Post to support the Nihon Kohden LifeScope G9 mobile monitor. The arm and device are mounted to a Mindray A7 Anesthesiology Workstation.



UL180EV7-W3-KUB-A1
A track mounted Ultra 180 arm supports the computer workstation on a Mindray A7 Anesthesiology Workstation.



Arm stows in 8" (20.3 cm)

ICW's versatile Ultra 180 arm configured with a Drop Post to support a Nihon Kohden LifeScope G9 mobile monitor.

- Vital patient monitoring services can be easily moved to bedside using ICW's new medical device mount attached to a Mindray A7 Anesthesiology Workstation.
- ICW's device mount can be adapted to fit a wide variety of medical and manufacturing/industrial devices.
- Devices are securely attached to the post with a specially-designed clamp.
- Devices can be adjusted up to 9" (23 cm) vertically on the Drop Post.
- The Ultra 180 Device Mount with Drop Post can support 32-40 lbs. (14.5-18 kg).

For more information about this new device mount, please contact your ICW representative, visit www.icwusa.com or call one of our product specialists at **1-800-558-4435**.

Kelly Tingley
General Manager
ICWUSA.com, Inc.

SPECIFICATIONS

Vertical adjustment	16" (40.6 cm)
Reach as shown	30" (76 cm)
Horizontal rotation at mount	180°
Horizontal rotation at A1 joint	up to 270°
Stowed depth	8" (20.3 cm)
Load capacity	32-40 lbs. (14.5-18 kg)